

SOP: Create a new EU standard test series by copying or cloning based on an existing EU test series. The EU test series clone can then be backdated in terms of validity, e.g. by 2, 3 or 4 years.

Step 1: Activate the extended test series management mode in the programme options.

2. call up the test series management

3. select the base EU test series e.g. 402

4. Copy the EU test series 402 to e.g. 502.

5. Select a historical date "valid from" in the date selection control e.g. 2 years before "today" or alternative 01.01.2020.

6. important: adjust test substances, vehicle and concentration if necessary

7. last important step: save the new EU test series 502.

8. the EU Test Series 502 can now also be used for data entry of patch test results for cases with a examination date from 01.01.2020 until today.